

ENT Endoscope Storage Cabinet

Stainless steel, 6 ENT endoscope capacity, wall mount storage cabinet with key lock



Contact us today!

+1.508.347.5100

sales@optim-llc.com

Optional LED-UV disinfection strip helps to keep scopes disinfected!



HEPA Filter Option

HEPA Filtered Positive Pressure with two high volume intake fans and one exhaust fan that activate automatically upon closing of doors. Change HEPA filters without opening cabinet, save time and ensure no cycle interruptions

Note: Installing the HEPA filter reduces scope storage capacity from 6 to 5 endoscopes.

Features & Benefits

Storage— Stores up to 6 ENT endoscopes

Vertical Design— Convenient access, while vertical position keeps optical insertion shaft in ideal position for protection and storage

Compatibility— Ideal for endoscopes without integral light guide cable or removable video/light processor cord

Ease of Use— Glass paneled doors for easy viewing of endoscope inventory, air flow vents on the side of the cabinet.

Quality— Industry best 12 year limited manufacturer warranty

Specifications

Cabinet dimensions: 30”H x 24”W x 12”D

Cabinet weight: 60 lbs

Contact us today to discuss your storage needs and discover the full range of solutions we offer.



64 Technology Park Road
Sturbridge, MA 01566 USA

+1.508.347.5100
sales@optim-llc.com

OPTIM'S PRODUCTS ARE
PROUDLY MANUFACTURED
IN THE USA