



## M-CAM HD USB Camera with TIMS Patient Management Software

*Designed to empower physicians with the visibility & mobility they need*

The add-on M-CAM HD USB Camera System makes recording video and taking still images easy and convenient. Simply couple to your flexible or rigid endoscope, plug the camera into the USB port on your PC or laptop and you are ready to go. It adds a new dimension of functionality to your endoscope, aiding in diagnosis by colleagues.

The easy to use, reliable and efficient TIMS Patient Management software includes patient database for visit comparison, image annotation, report creation for easy sharing and audio and video recording.



**Contact us today!**

+1.508.347.5100

[sales@optim-llc.com](mailto:sales@optim-llc.com)

### M-CAM HD USB Plug and Play Camera

**Easy to use**—USB Plug & Play, intuitive design, compatible with any qualified Windows-based computer, automatic white balance

**Compatibility**—Compatible with flexible and rigid endoscopes

**Image & Video Capture** — Easily capture and view high-quality images and videos using the Windows Camera app

**Zoom Coupler** — Provides precise image focusing with an adjustable 18 mm to 35 mm focal zoom for optimal image framing and clarity.

|                 |                      |
|-----------------|----------------------|
| Sensor          | CMOS Rolling Shutter |
| Shot Resolution | 2560 x 1920          |
| Weight          | 145g                 |
| Cable Length    | 3 Meters             |
| Protection Type | IP68 Waterproof      |

full HD



### TIMS Patient Management Software

**Compatibility**—Record studies at 30fps in full color at maximum resolution with synced audio

**Study Storage Options** —Archive studies to any PACS, VNA and EMR with the FULL TIMS Version or store directly to device with the TIMS LITE Version

**Productivity**— Quickly review video/image for annotation, reporting and more

**HIPPA Compliant**—Keep workflow digital. No memory sticks, CDs or DVDs to store and transport files

#### SIMPLE USB CONNECTION



Connect camera to laptop




Launch Windows Camera or TIMS software



Begin capturing images & video

Optim's ENTity XL Nasopharyngoscope paired with the M-CAM Camera connects directly to a laptop for clear, real-time image viewing.



64 Technology Park Road  
Sturbridge, MA 01566 USA

+1.508.347.5100  
[sales@optim-llc.com](mailto:sales@optim-llc.com)

Optim is ISO13485:2016 Certified

DOC-029, Rev A